

# nCounter™ GX Human GPCR Kit

Package insert for the nCounter™ GX Human GPCR Kit

## Product Description

The nCounter™ GX (Gene Expression) Human GPCR Kit is designed to produce an ultra-sensitive and highly multiplexed method for gene expression profiling. The nCounter Gene Expression Assay component consist of the CodeSet and all of the materials and reagents needed to perform the assay. Each of the sequences assayed by the CodeSet is Targeted by two different probes, the Reporter and the Capture probes.

## Reference Documents

- *nCounter™ Gene Expression Assay Manual - Total RNA and Cell Lysate Protocols*
- *nCounter™ CodeSet Design Manual*
- *nCounter™ Analysis System Brochure*
- *nCounter™ Prep Station User Manual*
- *nCounter™ Digital Analyzer User Manual*
- *nCounter™ Data Analysis Guidelines for Gene Expression*

## nCounter™ GX Human GPCR CodeSet

Component	Volume per Tube	Assays per Tube
Reporter ProbeSet (439 genes and controls)	130 µL	12
Capture ProbeSet	70 µL	12

Storage Conditions: Up to 6 months at -80°C.

The nCounter CodeSet kit also includes a USB drive that contains:

- Three read-only folders for Cartridge Definition files (CDFData), Reporter Library files (RLFData), and Reporter Code Count files (RCCData)
- A CDF template
- A PDF copy of the *nCounter™ Gene Expression Assay Manual - Total RNA and Cell Lysate Protocols*, the *nCounter™ Prep Station User Manual*, the *nCounter™ Digital Analyzer User Manual*, an RCC collector tool, *nCounter™ Data Analysis Guidelines for Gene Expression* and this product insert.
- Material Safety Data Sheet(s)

## Testing and Quality Control

All nCounter CodeSets are functionally tested.

## Specifications

Level of Multiplexing	439 Gene Targets per Assay <ul style="list-style-type: none"><li>• 431 GPCR gene targets</li><li>• 8 internal reference gene targets</li></ul>
Recommended Amount of Starting Material	100ng of total RNA, or lysate from ~10,000 cells
Sample Types Supported	Total RNA, Cell Lysates in GITC, FFPE derived total RNA and PAXgene lysed whole blood
Reaction Volume	30 µL
Limit of Detection	0.5fM spike-in control (~1 copy per cell); 90% of the time
Fold Change Sensitivity	>1.5 fold (>5 copies per cell) >2 fold change (>1 copy per cell)
Hybridization Spike Correlation	R <sup>2</sup> ≥ 0.95
Linear Dynamic Range	7 X 10 <sup>5</sup> Total Counts
nCounter Prep Station Throughput	12 samples <2.5 hours
nCounter Digital Analyzer Throughput	12 samples/4 hours (up to 72 samples per day unattended running in continuous mode)
Controls	6 Positive Controls 8 Negative Controls

## nCounter™ Master Kit Components

All reagents and consumables required to process nCounter Gene Expression Assays are supplied in the nCounter Master Kit. Kits for 12, 24, 48 assays are available. Orders for 96 Assays will be shipped two 48 assay Master Kits. Master Kit components include:

- nCounter Cartridges
- nCounter Prep Plates
- nCounter Prep Packs

The nCounter Master Kit components are shipped in separate boxes under different shipping conditions<sup>1</sup>.

### nCounter™ Cartridge

Twelve samples can be processed in a single nCounter Cartridge.

Component	Assay Pack	Cartridges per Pack
nCounter Cartridge	12	1
	24	2
	48	4

Storage Conditions: Store for up to 6 months at -20°C.

### nCounter™ Prep Plates

The Prep Plates are foil sealed 96-well plates that are used by the Prep Station. Do not remove the foil seals.

Component	Assay Pack	Plates per Pack
Contains wash buffers, magnetic beads, immobilization, and imaging reagents.	12	4
	24	4
	48	8

Storage Conditions: Store for up to 6 months at 4°C.

<sup>1</sup> Individual components may be purchased separately if required.

## nCounter™ Prep Pack

Component	Assay Pack	Units per Pack
Racked Tips & Foil Piercers	12	2 racks
	24	2 racks
	48	4 racks
12-Tube Strips	12	4 strips
	24	8 strips
	48	14 strips
Tip Sheaths	12	2 sheaths
	24	4 sheaths
	48	8 sheaths
12-Strip Tube Caps	12	4 strips
	24	8 strips
	48	14 strips
Cartridge Well Seals	12	2 seals
	24	3 seals
	48	6 seals
Hybridization Buffer*	12	1 tubes
	24	1 tubes
	48	2 tubes

\* If Hybridization buffer has been stored in cold conditions and precipitate is observed, warm tubes at 37°C until salts have dissolved.

Storage Conditions: Store at room temperature. No expiration date.

## Safety Information

Material Safety Data Sheets are included on the USB drive packaged with the nCounter CodeSets. Safety information for operation of the nCounter System is included in the *nCounter™ Prep Station User Manual* and the *nCounter™ Digital Analyzer User Manual*.

## Precautions

1. The NanoString nCounter System and Consumables are for Research Use Only; not for Diagnostic Procedures.
2. Avoid microbial contamination, which may negatively affect the quality of the results.
3. All biological specimens and materials should be handled as if the potential exists of transmitting infectious agents and disposed of with proper precautions in accordance with federal, state and local regulations. This includes adherence to the OSHA Blood Borne Pathogens Standard (29 CFR 1910.1030) for samples derived from blood and other sources governed by this act.
4. Never pipet by mouth.
5. Avoid specimen contact with skin and mucous membrane and always wear gloves.
6. Exercise caution when handling and disposing of carcinogenic reagents.
7. Avoid cross-contamination of samples which may negatively affect results quality.

## Contact Information

### NanoString Technologies®, Inc.

530 Fairview Ave N, Suite 2000

Seattle, WA 98109

Phone: 888-358-NANO (888-358-6266)

Fax: 206-378-6288

Email: [support@nanostring.com](mailto:support@nanostring.com)

[www.nanostring.com](http://www.nanostring.com)

## Limited License

Subject to the terms and conditions of the nCounter™ Analysis System contained in the product quotation, NanoString grants you a limited, non-exclusive, non-transferable, non-sublicensable, research use only license to use the proprietary nCounter Analysis System only in accordance with the manual and other written instructions provided by NanoString. Except as expressly set forth in the terms and conditions, no right or license, whether express, implied or statutory, is granted by NanoString under any intellectual property right owned by, or licensed to, NanoString by virtue of the supply of the proprietary nCounter Analysis System. Without limiting the foregoing, no right or license, whether express, implied or statutory, is granted by NanoString, to use the nCounter Analysis System with any third party product not supplied or licensed to you by NanoString, or recommended for use by NanoString in a manual or other written instruction provided by NanoString.

## Patents

The manufacture, use and/or sale of this product may be subject to one or more patents or pending patent applications owned by, or licensed to, NanoString Technologies, Inc. ("NanoString") including, but not limited to, one or more of the following patents: [US Patent Nos. 6,844,028 and 7,067,194 controlled by Accelr8 Technology Corporation, and the following patents: US Patent Nos. 5,696,157, 6,130,101, 6,974,873 and 6,977,305, and the following foreign patents: AU 0750380, AU 2001294859, EP 98948483.7, EP 0966458, EP 01975541.2, CA 2,272,403 and JP 11-519270, controlled by Molecular Probes, Inc.]

## Trademarks

NanoString Technologies, NanoString, nCounter and Molecules That Count are registered trademarks or trademarks of NanoString Technologies, Inc., in the United States and/or other countries. All other trademarks and/or service marks not owned by NanoString that appear in this package insert are the property of their respective owners.

## Copyright

© 2009 NanoString Technologies, Inc. All rights reserved.