

# #2441 NanoString 3D Biology™ Technology: Simultaneous digital counting of DNA, RNA, and Protein

Chris Lausted<sup>1</sup>, Chung-ying Huang<sup>2</sup>, Yong Zhou<sup>1</sup>, Gary Geiss<sup>2</sup>, Gavin Meredith<sup>2</sup>, Qian Mei<sup>2</sup>, Jinho Lee<sup>3</sup>, Gokhan Demirkan<sup>2</sup>, Dwayne Dunaway<sup>2</sup>, Nathan Elliott<sup>2</sup>, Dae Kim<sup>2</sup>, Elizabeth Manrao<sup>2</sup>, P. Martin Ross<sup>2</sup>, Warren Carter<sup>2</sup>, Lisa Bogatzki<sup>2</sup>, Sarah Warren<sup>2</sup>, Christopher P Vellano<sup>3</sup>, Karina A Eterovic<sup>3</sup>, Ping Song<sup>3</sup>, Lin-ya Tang<sup>3</sup>, Gloria Fawcett<sup>3</sup>, Tae-Beom Kim<sup>3</sup>, Ken Chen<sup>3</sup>, Gordon Mills<sup>3</sup>, Joseph Beechem<sup>2</sup>, and Leroy Hood<sup>1</sup>

<sup>1</sup> Institute for Systems Biology, Seattle WA. <sup>2</sup> NanoString Technologies, Seattle WA. <sup>3</sup> University of Texas M.D. Anderson Cancer Center, Houston TX.




530 Fairview Ave North Seattle, WA 98109

## Abstract

Development of improved cancer diagnostics and therapeutics requires detailed understanding of the genomic, transcriptomic, and proteomic profiles in the tumor microenvironment. Current technologies can excel at measuring a single analyte-type, but it remains challenging to simultaneously collect high-throughput DNA, RNA, and protein data from small samples.

We have developed an approach that uses optical barcodes to simultaneously profile DNA, RNA, and protein from as little as 5 ng DNA, 25 ng RNA, and 250 ng protein or just two 5 μm FFPE slides, and simplifies data analysis by generating digital counts for each analyte.

## 3D Biology™ Technology Overview




**FIGURE 1.** Principles of the multiplex 3D Biology assay. Up to 800 targets from any combination of DNA, RNA, and protein can be interrogated from a single sample by labeling with optical barcodes using either direct hybridization of barcode-conjugated oligonucleotides (for DNA and RNA) or via oligo-conjugated primary antibodies (for protein). Analyte types are then pooled into a single reaction and imaged on the nCounter® Analysis system. Quantification of the barcodes generates a representative count of the number of molecules of target that was present in the original sample.

## 3D Biology™ Technology Highlights

- Low sample input**  
5 ng DNA  
25 ng RNA  
250 ng Protein  
As few as 5,000 cells or two 5 μm FFPE sections
- 3 analytes measured from the same biological sample**  
104 SNVs from 25 key driver genes  
Up to 770 mRNA targets  
30+ Proteins
- One instrument provides digital readout of SNV, gene expression, and protein profiling from minimal sample, including FFPE tumor tissue**

## DNA: Single Nucleotide Variant (SNV) Detection




A single nucleotide mismatch is thermodynamically unstable and disrupts probe hybridization.

Hybridization complexes are captured on the imaging surface and digitally quantified to reveal relative differences between wild type (WT) and variant alleles.

**2B**


Probe Name	Codon 12	Codon 13	COSM16	COSM20	COSM21	COSM22	COSM27	COSM27	COSM27	COSM32	COSM17	COSM12	COSM18	COSM29	COSM28
Ref (exon 2)	G	G	T	G	G	C	A	A	A	A	A	A	A	A	A
COSM16 (G12C)	T		14.5%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%	0.2%
COSM20 (G12V)	T		0.2%	36.9%	0.7%	0.4%	0.5%	0.5%	0.5%	0.5%	0.7%	0.5%	0.5%	0.5%	0.5%
COSM21 (G12D)	A		0.4%	0.4%	18.5%	0.3%	0.4%	0.4%	0.3%	0.4%	0.5%	0.3%	0.3%	0.3%	0.3%
COSM22 (G12A)	C		0.2%	0.3%	0.2%	29.5%	0.2%	0.2%	0.2%	0.3%	0.2%	0.2%	0.2%	0.2%	0.2%
COSM27 (G13C)		T				17.2%	0.1%	0.1%	0.1%	0.1%	0.1%	0.1%	0.1%	0.1%	0.1%
COSM32 (G13D)			A			0.3%	0.3%	0.4%	0.4%	0.3%	0.4%	0.4%	0.4%	0.3%	0.4%
COSM17 (G12S)	A		0.6%	0.6%	0.6%	0.4%	0.5%	0.5%	0.5%	0.6%	0.4%	0.5%	0.4%	0.5%	
COSM12 (G12F)	T	T	0.2%	0.2%	0.1%	0.1%	0.2%	0.2%	0.1%	0.2%	0.1%	0.1%	0.1%	0.1%	
COSM18 (G12R)	C		0.2%	0.3%	0.3%	0.2%	0.2%	0.2%	0.2%	0.3%	0.2%	0.2%	0.2%	0.2%	
COSM33 (G13A)			C			0.4%	0.4%	0.3%	0.3%	0.4%	0.3%	0.5%	0.3%	0.4%	
COSM29 (G13R)			C			0.2%	0.2%	0.1%	0.1%	0.1%	0.1%	0.2%	0.1%	0.1%	
COSM28 (G13S)			A			0.3%	0.3%	0.3%	0.3%	0.3%	0.3%	0.3%	0.3%	0.3%	

**FIGURE 2. A)** Diagram of SNV probe architecture to differentiate WT from mutant allele. **B)** KRAS exon 2, codons 12 and 13, contain several known SNV. gDNA standards containing several of these mutations (COSMIC ID) were pooled and tested using the nCounter® Vantage 3D™ SNV Assay. Specificity is determined by the percentage of counts for a given probe to its intended target compared to counts for all KRAS targets.




**FIGURE 3.** gDNA extracted from 45 different FFPE tumor samples with known SNV content (determined by NGS) was tested by nCounter Vantage 3D SNV Assay using a panel of 104 SNV probes. **True positive:** variant detected by NGS and nCounter. **Putative False Positive:** variant detected by nCounter but not NGS. **False Negative:** variant detected by NGS but not nCounter.


## Multiplexed RNA and Protein Analysis



As few as 20,000 cultured cells (50,000 primary cells) can be collected by the Universal Cell Capture Beads, and processed for both RNA and protein quantification.




mRNA and protein expression of Jurkat-specific CD3E and CD28 are highly correlated in titration of Jurkat into Karpas cells.



nCounter protein count linearity is higher than SRM mass spectrometry for CD3E and CD28 in Jurkat cell model.

## 3D Biology Application



**FIGURE 5.** Application of 3D Biology Technology. **(A)** BRAF<sup>WT</sup> and BRAF<sup>V600E</sup> cell lines, treated with DMSO or Vemurafinib, are profiled for 104 SNV, 30 protein, and 192 mRNA targets. Only data for genes with >100 digital counts shown. **(B)** Volcano plot of 8hour Vemurafinib treatment vs. DMSO, showing significant protein and mRNA fold-changes as a result of treatment.

## Methods

**SNV Assay:** Genomic DNA was purified from single, 5 μm FFPE tissue sections or 25,000 cultured cells using the Qiagen AllPrep kit. Purified DNA was quantified by Qubit fluorometric assay. Exactly 5 ng sample DNA was used as input in the SNV amplification reaction. DNA integrity scores were determined by Agilent TapeStation analysis.

**Cell collection and RNA/protein sample lysis:** PBMC samples (Centrifugation method): 500,000 PBMCs were collected by centrifugation at 400 RCF for 5 minutes. Cell pellets were lysed in 77 μL Qiagen buffer RLT (8500 cells/μL RLT), 1.5 μL, or 10,000 cell equivalents, were input into the nCounter hybridization. Universal Cell Capture method: 50,000 PBMC were collected using the nCounter Vantage 3D Universal Cell Capture Kit. Briefly, PBMC cell suspension is incubated with Universal Cell Capture beads (anti-β2M conjugated magnetic beads) and collected on a 96-well plate magnet, supernatant is removed and Fc Receptor Blocking solution added. An aliquot of 50,000 blocked cells is collected by magnet and a portion are incubated with Antibody mix, collected and lysed in buffer LH with another portion dedicated to direct lysis with buffer LH for RNA analysis. 10,000 cell equivalents of lysate are added to the RNA hybridization, while 800 cell equivalents are added for protein hybridization.

**Cell line samples:** Jurkat and Karpas cell lines were cultured separately, harvested, and mixed at indicated percentages to create mock fold-change samples. BRAF<sup>WT</sup>, BRAF<sup>V600E</sup> heterozygous, and BRAF<sup>V600E</sup> homozygous mutant cell lines (SKMEL2, SKMEL5, SKMEL28) were cultured from 10min to 24 hours in presence of 0.1% DMSO or 250 nM Vemurafinib (1mM stock in DMSO, Selleckchem Houston, TX). NanoString has validated the RNA:protein assay with as little as 20,000 cultured cells collected by Universal Cell Capture.

## Conclusions

NanoString 3D Biology Technology enables highly sensitive and reproducible digital profiling of DNA, RNA, and protein from small amounts of cells and tumor FFPE tissue. We have presented data demonstrating several key attributes of this experimental workflow:

- ✓ High SNV probe specificity that performs comparably to the NGS benchmark
- ✓ Protein detection that correlates highly with both flow cytometry and SRM mass spectrometry
- ✓ Recapitulation of known genotype-dependent gene expression and protein phosphorylation changes in a melanoma cancer model

