


PBMC Tips and Tricks

WHITE PAPER

Emer Clarke, Ph.D.
Gary dos Santos



PBMC Tips and Tricks: How to make the most of this precious sample type

Emer Clarke, Ph.D., and Gary dos Santos

Introduction

In recent years, there has been a movement from the use of cell lines to primary cells to address a number of biological questions including, but not restricted to, target validation, toxicity, immunology and cytokine storm risk. Although cell lines are easy to use, predictable in their growth patterns, and amenable to high-throughput assays, they are generally limited in simple end points such as cell death or cell proliferation. Cell lines usually represent a mature cell type from only one lineage and are prone to genetic instability. Therefore, the translation of information from these assays may not correlate well with the *in vivo* situation. As our awareness increases as to the many players in the immune system, so too does our understanding of the function of these cells. For example, some cancer cells are known to “exhaust” T cell activity¹, thus dampening the immune response, whereas other engineered T cells are being evaluated to overcome their loss of activity in the tumor microenvironment (Figure 1).^{2,3}




FIGURE 1 Example of T cell inactivation through PD1/PD-L1

As the scientific community addresses more specific queries about individual cell function, both alone and in the context of other cells, the need for more representative cellular platforms is required. Peripheral blood mononuclear cells (PBMCs) provide one such platform. The ability to acquire these cells has become easier, with many reputable vendors throughout the United States. The added advantage of accessing blood from both normal donors and patients with various diseases opens the possibility of researching the similarities or differences in target number, signaling pathways, or molecular response to drug candidates in these populations. The ability to “bank” many identical vials of cells from a given donor allows for comparative analyses of test compounds over time, whereas the ability to access PBMCs from multiple donors can assess if a test compound could have varying molecular effects on a population. However, working with PBMCs can be more difficult than cell lines. Part of this white paper hopes to present several pointers to aid the generation of viable and functional starting cell populations.

What are PBMCs?

PBMCs are the mononuclear cells within the blood. This includes: lymphocytes (T cells, B cells, Natural Killer (NK) cells, and monocytes).⁴ The relative distribution of these cell populations may be influenced by age, ethnicity, or gender.⁵ From whole blood, which additionally contains red blood cells, platelets, and granulocytes, PBMCs can be prepared by layering the cells on a density gradient (ficoll) and centrifuging the tube (Figure 2). Denser cells will sink to the bottom, whereas the PBMCs can be found in a layer between the ficoll and the plasma.⁶ PBMCs may also be acquired from apheresis, a process in which blood from a given donor or patient is passed through an apparatus that separates out one particular constituent (PBMCs in this case) and returns the remainder of the blood components to the donor.^{7,8,9}




FIGURE 2 Representation of PBMC preparation before and after density gradient centrifugation.

Contaminating Cells in PBMC Preparations

Despite best efforts by density centrifugation or apheresis, there can be contaminating cells, some of which cause more problems than others. Contaminating red blood cells are typically not too problematic. If the PBMCs are cryopreserved and stored in the gaseous phase of liquid nitrogen prior to use, the red blood cells do not survive these temperatures and are no longer a problem upon thawing of the sample. If fresh PBMCs are required in a given assay, contaminating red blood cells would result in an erroneous cell count. This can be avoided by the use of a calibrated cell counter or by dilution of the PBMC samples in glacial acetic acid (which lyses red blood cells) for a manual cell count.

Granulocyte contamination causes more challenges. Granulocytes have the shortest half-life of all blood cells, so as blood ages (between collection time and use), they begin to change density, which may reduce the effectiveness of a density gradient separation. In addition, with a half-life of approximately 9 hours, they can die and release products into the PBMCs, some of which are known to alter the function of certain populations including T-cells.¹⁰ Granulocytes also do not survive cryopreserving and storage in the gaseous phase of liquid nitrogen, so upon thawing, a PBMC population contaminated with granulocytes may yield a poor cell viability.

Use of Fresh vs. Cryopreserved Cells

For many applications, cryopreserved PBMCs are an ideal cell population. The ability to have a cell stock ready for use has obvious advantages. In order to ensure intended assays work optimally, it is important to thaw cryopreserved PBMCs rapidly.

Once removed from the gaseous phase of liquid nitrogen, the vial(s) of cells should be placed in a 37°C water bath, until they are almost completely thawed, then transferred to a basal medium containing 5-10% FBS (or a medium of choice prepared for the assay) and washed by centrifugation at approximately 1200 rpm for 10 minutes. Following decanting of the supernatant and prior to the resuspension of the cells, 50-100 µL DNase may be added. This aids in preventing cell clumping but should not be used for any assays measuring DNA or that utilize DNA for assay readout. Following resuspension of cells in the medium of choice, a cell and viability count should be performed. Once the PBMCs have been diluted to the appropriate cell concentration, the assays should be set up immediately. Samples with a low viability (i.e. many dead cells) are unlikely to yield optimal assay results.

Where use of fresh PBMCs is required, the fresh blood needs to be processed as rapidly as possible. Typically, the blood is diluted in basal medium, layered on ficoll and the PBMC population recovered following density gradient separation. Following washing of the PBMCs, DNase may be added again to prevent any clumping of cells, but again should not be used in assays that require DNA analysis or DNA for assay readout. A cell and viability count should be performed. Since dead cells have a different density than their live counterparts, the ficoll separation has the advantages of removing the dead cells in the process so the viability of PBMCs is typically > 95% when prepared freshly.

DO'S AND DON'Ts WITH PBMCs:

DO...

Do use all the appropriate safety equipment (laminar flow hoods and gloves) when processing or working with human blood or PBMCs.

Do process fresh samples as soon as possible after receipt.

Do perform a viability count, especially on thawed PBMC samples.

Do have all other reagents for the culture ready to go so the PBMCs (fresh or thawed) are used immediately upon preparation.

DON'T...

Don't leave cryopreserved cells in a water bath for extended time periods when thawing (the vial should be removed when there is still evidence of some ice crystal).

Don't use DNase if you are planning on using the cells in a DNA-based assay such as the nCounter® platform for protein analysis.

Don't store cryopreserved cells in freezers (-20°C, -80°C) as the functionality of some cell populations is compromised.

Assay formats may vary depending on the scientific question, the desired readouts and the cell numbers required for the various readouts, but a 96-well plate format offers an excellent platform for many types of experiments with PBMCs.

Conclusion

In summary, PBMCs are readily available and provide a relevant cell population for addressing many biological queries. With a little additional attention to detail in cell preparation and cell storage, the data generated from viable PBMCs may provide a wealth of information about normal versus diseased and treated versus untreated patients, and information regarding potential drug candidates with an increased understanding of underlying mechanisms.

Author Profiles



Emer Clarke, Ph.D. is co-founder and Chief Science Officer at ReachBio Research Labs. She is instrumental in managing the scientific operations of the company including research and development, and client services. Previously, she was a member of the senior management team at StemCell Technologies Inc. as Scientific

Director of Contract Services and Manager of Education. Dr. Clarke earned her Ph.D. in hematology from Trinity College Dublin.



Gary dos Santos is the Director of Laboratory Operations at ReachBio Research Labs. With over 17 years in the biotech field, his expertise has been essential in helping expand the company's drug development platforms. Prior to becoming ReachBio's first employee, nine years ago, he worked as a scientist for

StemCell Technologies Inc., and Veterinary Pathology Laboratory Inc. He earned both his B.Sc. in biochemistry, and his M.Sc. in microbiology from McGill University in Montreal, Canada.



References

1. Sugimoto, T. & Watanabe T. (2016). Follicular Lymphoma: The Role of the Tumor Microenvironment in Prognosis. *J Clin Exp Hematop* 56(1):1-19.
2. Davila M.L. & Sadelain M. Biology and Clinical Application of CAR T Cells for B cell malignancies. *Int. J. Hematol* 104(1): 6-17.
3. Cherkassy L. et al (2016). Human CAR T Cells with cell-intrinsic PD-1 checkpoints blockade resist tumor-mediated inhibition. *J. Clin Invest.* 126 (8): 3130-3144.
4. Goff L. et al (1985). Normal values for the different classes of cenous blood mononuclear cells defined by monoclonal antibodies. *J Clin Path* 38(1): 54-59.
5. Tollerud D.J. et al (1989). The influence of age, race, and gender on peripheral blood mononuclear-cell subsets in healthy non-smokers. *J Clin Immunol* 9(3): 214-222
6. Casale T.B. & KAliner M. (1982). A rapid method for isolation of human mononuclear cells free of significant platelet contamination. *J Immunol Methods* 55(3): 347-353.
7. Mookerjee B.K. (1076). Influence of separation techniques on the distribution and function of lymphocyte population. A comparison on three techniques. *Transplantation* 22 (2); 101-107.
8. Snyder E.L. et al (2001). Ex vivo evaluation of PBMNCs collected with a new separator. *Transfusion* 41(7): 940-949.
9. Steininger P.A. et al (2014). First comparative evaluation of a new leukapheresis technology in non-cytokine stimulated donors. *Vox Sang* 106(3): 248-255.
10. McKenna K.C. et al (2009). Delayed processing of blood increases the frequency of activated CD11b+ CD15+ granulocytes which inhibit T cell function. *J Immunol Methods* (341 (1-2): 68-75.

NanoString Technologies, Inc.

530 Fairview Ave. North
Seattle, Washington 98109 USA

LEARN MORE

Visit www.nanostring.com/products/vantage-rna-protein to learn more about nCounter Vantage 3D™ Immune Cell Assays.

Tel: (888) 358-6266 | Fax: (206) 378-6288
www.nanostring.com | info@nanostring.com

SALES CONTACTS

United States: us.sales@nanostring.com
EMEA: europe.sales@nanostring.com
Asia Pacific & Japan: apac.sales@nanostring.com
Other Regions: info@nanostring.com